

TVRDOKOVOVÉ A DIAMANTOVÉ VRTÁČIKY PRE NSK TURBÍNKU S-MAX PICO



Diamantové vrtáčky – NSK

Model	PC1	PC2	PC3	PC4	PC5	PC6
bur						
REF (Packof 3 pc)	Y1001817	Y1001818	Y1001819	Y1001820	Y1001821	Y1001822
length total	16,8mm	16,8mm	17,5mm	16,0mm	14,0mm	15,5mm
length of diamond coating	1,9mm	0,4mm	9,0mm	7,0mm	4,0mm	6,0mm
diameter diamond coating	0,9mm	0,6mm	1,8mm	1,3mm	1,6mm	1,1mm
grit	Standard	Fein	Grob	Standard		
ISO code	237/009	001/006	198/018	165/013	170/016	288/011

Cena: 1 krabička(3 ks) 33 €



Tvrdokovové vrtáčky – Komet

H7



	5	5	5	5	5	5	5
Größe · Size	Ø 1/10 mm 006	008	009	010	012	014	018
L mm	1,2	1,6	1,7	1,7	1,8	1,8	2,4
US No.	329	330	-	331	332	333	-

FG kurz · FG short



500 313 232001 ...
H7.313. ...



008

max. 100 000 min⁻¹

S vysokou účinnosťou rezania a pre exkaváciu


Cena: 3,74 €/ks*

H59



	5
Größe · Size	Ø 1/10 mm 010
L mm	2,5

FG kurz · FG short



H59.313. ...



010

Pre otvorenie fisúr, minimálne invazívne

Cena: 5,83 €/ks*

Diamantové vrtáčky - Komet

8889 M
889 M

i

		5
Größe · Size	Ø 1/10 mm	007
L mm		2,7

8838 M
838 M

i

		5
Größe · Size	Ø 1/10 mm	007
L mm		2,7

8830 M
830 M

i

		5
Größe · Size	Ø 1/10 mm	012
L mm		2,7

FG kurz · FG short

889 M.313. ... **007**

Cena: 5,50 €/ks*

FG kurz · FG short

838 M.313. ... **007**

Cena: 5,50 €/ks*

FG kurz · FG short

830 M.313. ... **012**

max. 160 000 min⁻¹

Cena: 5,50 €/ks*

8953 AM
953 AM

i

		5
Größe · Size	Ø 1/10 mm	014
L mm		2,5

Cena: 5,50 €/ks*

8953 M
953 M

i

		5
Größe · Size	Ø 1/10 mm	014
L mm		2,0

Cena: 5,50 €/ks*

8830 RM
830 RM

i

		5
Größe · Size	Ø 1/10 mm	009
L mm		2,7

Cena: 5,50 €/ks*

FG kurz · FG short

953 AM.313. ... **014**

max. 160 000 min⁻¹

FG kurz · FG short

953 M.313. ... **014**

max. 160 000 min⁻¹

FG kurz · FG short

830 RM.313. ... **009**

max. 300 000 min⁻¹

801 UF
 801 EF
 8801
 801
 6801
 5801

		5	5	5	5	5	5	5	5	5	5	5	5	5
Größe · Size	Ø 1/10 mm	006	007	008	009	010	012	014	016	018	023	025	029	035

FG kurz · short


801.313. ...

806 313 001524 ...

010 014 018

Cena: 2,97 €/ks*

- **830 LEF**
- **8830 L**
- 830 L**
- **6830 L**
- **5830 L**



		5	5	5	5	5
Größe · Size	Ø 1/10 mm	010	012	014	016	018
L mm		4,0	4,0	5,0	5,0	5,0

FG kurz · short



806 313 234524 ...
830L.313. ...

				014		

Cena: 5,50 €/ks*

- **830 EF**
- **8830**
- 830**
- **6830**



		5	5	5	5	5	5
Größe · Size	Ø 1/10 mm	007	008	009	010	012	014
L mm		2,7	2,7	2,7	2,7	2,7	2,7

FG kurz · short



806 313 233524 ...
830.313. ...

					012		

Cena: 5,50 €/ks*

- **8835**
- 835**
- **6835**



		5	5	5	5	5	5
Größe · Size	Ø 1/10 mm	007	008	009	010	012	014
L mm		3,0	3,0	3,0	4,0	4,0	4,0

FG kurz · short




806 313 107524 ...
835.313. ...

					010	012	
						012	

■ = max. 100 000 min⁻¹

Cena: 5,50 €/ks*

- **858 UF**
- **858 EF**
- **8858**
- 858**
- **6858**



		5	5	5
Größe · Size	Ø 1/10 mm	010	014	016
L mm		8,0	8,0	8,0
Winkel · Angle α		2°	3,4°	3,9°

FG kurz · FG short




806 313 165524 ...
858.313. ...

					014	

◆ = max. 300 000 min⁻¹

Cena: 2,97 €/ks*

- **862 UF**
- **862 EF**
- **8862**
- 862**
- **6862**
- **5862**



		5	5	5	5
Größe · Size	Ø 1/10 mm	010	012	014	016
L mm		8,0	8,0	8,0	8,0

FG kurz · short



806 313 249534 ...
6862.313. ...


				012		

■ = max. 100 000 min⁻¹

◆ = max. 300 000 min⁻¹

Cena: 5,50 €/ks*

- **8877**
- 877**
- **6877**



		5	5	5	5	5
Größe · Size	Ø 1/10 mm	009	010	012	014	016
L mm		6,0	6,0	6,0	6,0	6,0

FG kurz · short



806 313 288514 ...
8877.313. ...

				012		
				012		

◆ = max. 160 000 min⁻¹

Cena: 3,58 €/ks*



	5	5	5	5	
Größe - Size	Ø 1/10 mm	010	012	014	016
L mm	8,0	8,0	8,0	8,0	

FG kurz - short

806 313 289514 ...	-	◆012	-	-
806 313 289524 ...	-	◆012	-	-
806 313 289534 ...	-	-	014	-

◇ = ⌀ max. 160 000 min⁻¹
◆ = ⌀ max. 300 000 min⁻¹



	5	5	5	5	
Größe - Size	Ø 1/10 mm	018	023	027	031
L mm	2,3	2,6	2,7	3,0	

FG kurz - short

8905.313. ...	◆018	◆023	◆027	◆031
905.313. ...	◆018	◆023	◆027	◆031

zum okklusalen Gestalten
for occlusal shaping

◆018 = ⌀ max. 300 000 min⁻¹
◆023 = ⌀ max. 160 000 min⁻¹
◆027-031 = ⌀ max. 140 000 min⁻¹

Pre okluzálne tvarovanie

Cena: 7,26 €/ks*

Cena: 5,50 €/ks*

*Poznámka: všetky vrtáčky sa objednávajú v min. množstve 5 ks (blister)

Diamantové vrtáčky – NTI (Germany)



801 (001)
Rund
Round
Rond



Größe/Size/taille ø 1/10 mm	007	008	009	010	012	014	016	018	021	023	025	029	035	042	
Länge - Length in mm - longueur															
Shank	ISO	REF													
FGM 806 313 001 524...	801- ... M- FGM			012			014			018					
FGM 806 313 001 534...	801- ... C- FGM			023											

Cena: 40 €/krabičku (5 ks)

805 (010)
Umgekehrter Kegel
Inverted Cone
Cône renversé



Größe/Size/taille ø 1/10 mm	008	009	010	012	014	016	018	021	023
Länge - Length in mm - longueur	1,0	1,0	1,0	1,2	1,6	1,6	1,6	1,9	2,5
Shank	ISO	REF							
FGM 806 313 010 524...	805- ... M- FGM			012		014			

Cena: 40 €/krabičku (5 ks)

806 (019)
Umgekehrter Kegel, mit Ansatz
Inverted Cone with Collar
Cône renversé, avec extension



Größe/Size/taille ø 1/10 mm	008	009	010	012	014	016	018	023
Länge - Length in mm - longueur	2,0	2,0	2,0	3,0	3,0	3,0	3,0	4,0
Shank	ISO	REF						
FGM 806 313 019 524...	806- ... M- FGM			010				

Cena: 40 €/krabičku (5 ks)

392 (465)
Interdental
Interproximal
Interdentaire



Größe/Size/taille ø 1/10 mm	016
Länge - Length in mm - longueur	5,0
Shank	ISO REF
FGM 806 313 465 524...	392 - ... M- FGM 016

**835
(107)** 

Zylinder
Flat End Cylinder
Cylindre à bout plat




Größe/Size/taille ø 1/10 mm	007	008	009	010	012	014	016	018
Länge-Length in mm - longueur	3,0	3,0	3,0	4,0	4,0	4,0	4,0	4,0
Shank <small>ISO</small>	<small>REF</small>							
FGM 806 313 107 524...	835- ... M- FGM	008	009	010	012	014		
FGM 806 313 107 534...	835- ... C- FGM	010						

Cena: 40 €/krabičku (5 ks)


**836
(110)** 

Zylinder
Flat End Cylinder
Cylindre à bout plat




Größe/Size/taille ø 1/10 mm	008	010	012	014	016	018	027	
Länge-Length in mm - longueur	6,0	6,0	6,0	6,0	6,0	6,0	6,0	
Shank <small>ISO</small>	<small>REF</small>							
FGM 806 313 110 544...	836- ... SC-FGM	014						

Cena: 40 €/krabičku (5 ks)

**837
(111)** 

Zylinder
Flat End Cylinder
Cylindre à bout plat



Größe/Size/taille ø 1/10 mm	010	012	014	016	018
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0
Shank <small>ISO</small>	<small>REF</small>				
FGM 806 313 111 534...	837- ... C- FGM	014			

Cena: 40 €/krabičku (5 ks)

**845
(170)** 


Konus
Flat End Taper
Cône à bout plat



Größe/Size/taille ø 1/10 mm	009	010	012	014	016
Länge-Length in mm - longueur	3,0	4,0	4,0	4,0	4,0
Shank <small>ISO</small>	<small>REF</small>				
FGM 806 313 170 524...	845- ... M- FGM	010	012		

**847
(172)** 

Konus, flach
Flat End Taper
Cône à bout plat




Größe/Size/taille ø 1/10 mm	010	012	014	016	018	023
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0
Shank <small>ISO</small>	<small>REF</small>					
FGM 806 313 172 544...	847- ... SC-FGM	016				
FGM 806 313 172 534...	847- ... C- FGM	016				


Cena: 40 €/krabičku (5 ks)

**847KR
(546)** 


Konus, Kante rund
KR Taper, Modified Shoulder
Bout arrondi



Größe/Size/taille ø 1/10 mm	016	018	023
Länge-Length in mm - longueur	8,0	8,0	8,0
Shank <small>ISO</small>	<small>REF</small>		
FGM 806 313 546 544...	847KR- ...SC-FGM	016	

**848
(173)** 

Konus, flach
Flat End Taper
Cône bout plat




Größe/Size/taille ø 1/10 mm	010	012	014	016	018	023
Länge-Length in mm - longueur	10,0	10,0	10,0	10,0	10,0	10,0
Shank <small>ISO</small>	<small>REF</small>					
FGM 806 313 173 524...	848- ... M- FGM	014				

Cena: 40 €/krabičku (5 ks)

**850
(199)**

Konus, rund
Round End Taper
Cône rond



Größe/Size/taille ø 1/10 mm	010	012	014	016	018	023	031
Länge-Length in mm - longueur	10,0	10,0	10,0	10,0	10,0	10,0	10,0
Shank ISO REF							
FGM 806 313 199 534... 850- ...C- FGM							014

Cena: 40 €/krabičku (5 ks)

**850KR
(553)**


Konus, Kante rund
KR Taper, Modified Shoulder
Bout arrondi



Größe/Size/taille ø 1/10 mm	014	018
Länge-Length in mm - longueur	10,0	10,0
Shank ISO REF		
FGM 806 313 553 524... 850KR-...M- FGM	018	

**855
(197)**

Konus, rund
Round End Taper
Cône rond




Größe/Size/taille ø 1/10 mm	012	014	016	018	025
Länge-Length in mm - longueur	6,0	6,0	6,0	6,0	7,0
Shank ISO REF					
FGM 806 313 197 534... 855- ...C- FGM					012

Cena: 40 €/krabičku (5 ks)

**858
(165)**

Konus, spitz
Needle
Cône pointu



Größe/Size/taille ø 1/10 mm	010	012	014	016	018
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0
Shank ISO REF					
FGM 806 313 165 524... 858- ...M- FGM					014
FGM 806 313 165 534... 858- ...C- FGM					014
FGM 806 313 165 504... 858- ...SF- FGM					014

**856
(198)**

Konus, rund
Round End Taper
Cône rond




Größe/Size/taille ø 1/10 mm	009	010	012	014	016	018	021	025	033
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0	8,0	8,0	9,0
Shank ISO REF									
FGM 806 313 198 524... 856- ...M- FGM									012
FGM 806 313 198 544... 856- ...SC-FGM	012	014	016	018	021	025			
FGM 806 313 198 534... 856- ...C- FGM	012	014	016	018	021	025			

Cena: 40 €/krabičku (5 ks)

**859
(166)**

Konus, spitz
Needle
Cône pointu



Größe/Size/taille ø 1/10 mm	010	012	014	016	018	021	
Länge-Length in mm - longueur	10,0	10,0	10,0	10,0	10,0	10,0	
Shank ISO REF							
FGM 806 313 166 524... 859- ...M- FGM							016

Cena: 40 €/krabičku (5 ks)

**860
(247)**


Flamme
Flame
Flamme



Größe/Size/taille ø 1/10 mm	010	012	014	016
Länge-Length in mm - longueur	4,0	5,0	5,0	5,0
Shank ISO REF				
FGM 806 313 247 534... 860- ...C- FGM				012

**862
(249)**

Flamme
Flame
Flamme



Größe/Size/taille ø 1/10 mm	009	010	012	014	016	018	021
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0	8,0
Shank ISO REF							
FGM 806 313 249 524... 862- ...M- FGM							016
FGM 806 313 249 534... 862- ...C- FGM		010	012	014			

Cena: 40 €/krabičku (5 ks)

**863
(250)**

Flamme
Flame
Flamme



Größe/Size/taille ø 1/10 mm	010	012	014	016	018	021
Länge-Length in mm - longueur	10,0	10,0	10,0	10,0	10,0	10,0
Shank ISO REF						
FGM 806 313 250 534...	863- ... C- FGM		012			

Cena: 40 €/krabičku (5 ks)

**878
(289)**

Torpedo
Modified Beveled Cylinder
Torpille




Größe/Size/taille ø 1/10 mm	008	009	010	012	014	016	018
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0	8,0
Shank ISO REF							
FGM 806 313 289 534...	878- ... C- FGM		012		014		

Cena: 40 €/krabičku (5 ks)

**879K
(299)**

Torpedo, konisch
Curettage
Torpille conique, Curettage




Größe/Size/taille ø 1/10 mm	012	014	016	018	021	023
Länge-Length in mm - longueur	10,0	10,0	10,0	10,0	10,0	10,0
Shank ISO REF						
FGM 806 313 299 524...	879K- ... M- FGM		016			

Cena: 40 €/krabičku (5 ks)

**878K
(298)**


Torpedo, konisch
Curettage
Torpille conique, Curettage



Größe/Size/taille ø 1/10 mm	012	014	016	018	021	023
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0
Shank ISO REF						
FGM 806 313 298 534...	878K- ... C- FGM		016			


**885
(130)**

Zylinder, spitz
Beveled Cylinder
Cylindre pointu



Größe/Size/taille ø 1/10 mm	008	010	012	014	016	018
Länge-Length in mm - longueur	8,0	8,0	8,0	8,0	8,0	8,0
Shank ISO REF						
FGM 806 313 130 534...	885- ... C- FGM		012			


Cena: 40 €/krabičku (5 ks)



**Kiddies FGM
REF SET-1669**

Cena setu: 85 €

Körnung, Grit, Grain
ISO 514
62 - 74 µm



Größe/Size/taille ø 1/10 mm	009	007	008	010	008	007	007	012	007	007	007	007
Länge-Length in mm - longueur	1,0	-	1,0	1,0	3,0	2,0	3,5	6,0	2,0	2,0	2,0	3,5
Shank ISO REF												
FGM 806 313 277 514...	379- ... M- KID		009									
FGM 806 313 001 514...	801- ... M- KID		007									
FGM 806 313 012 514...	805- ... M- KID		008									
FGM 806 313 010 514...	805- ... M- KID		010									
FGM 806 313 108 514...	835- ... M- KID		008									
FGM 806 313 137 514...	838- ... M- KID		007									
FGM 806 313 138 514...	838L- ... M- KID		007									
FGM 806 313 171 514...	846- ... M- KID		012									
FGM 806 313 161 514...	852- ... M- KID		007									
FGM 806 313 246 514...	860- ... M- KID		007									
FGM 806 313 495 514...	883- ... M- KID		007									
FGM 806 313 496 514...	889- ... M- KID		007									

opt. 160.000 rpm

Uvedené ceny sú s 20% DPH.