

**NSK**

RECHARGEABLE MICROMOTOR

*Viva-mate 3*

high performance & compact



Powerful Partners®

# Viva-mate 3

## Rechargeable Micromotor

Viva-mate 3 is a rechargeable and portable coreless micromotor system with exceptional power. Including integrated battery, the entire control unit is compact and weighs only 600g which offers great portability. Quick one hour battery charge allows one hour of continuous use. Attaching different handpieces, Viva-mate 3 is capable for various use such as polishing, general cutting, endodontic treatment and prosthetic adjustment procedures.

### Features

- High performance micromotor
- Quick one hour charge
- Charging indicator
- Foot control (Optional)

### Specifications

- Power supply : Ni-Cd battery
- Charging time : 1 hour
- Speed : 2,000~30,000 min<sup>-1</sup>
- Unit weight : 600 g
- Unit dimensions : W72 x D163 x H39 (mm)



	MODEL	ORDER CODE
	EX-6B	H260-001
	<ul style="list-style-type: none"> <li>• Included</li> <li>• 1:1 Direct Drive • For HP burs (ø2.35)</li> </ul>	
	AR-ER4 (S)	Y110-140
	<ul style="list-style-type: none"> <li>• Optional</li> <li>• 4:1 Reduction • For Snap-on Cups</li> </ul>	
	TEP-ER10	Y110-145
	<ul style="list-style-type: none"> <li>• Optional</li> <li>• 10:1 Reduction • For Hand Files</li> </ul>	
	Ti95EX	C487
	<ul style="list-style-type: none"> <li>• Optional</li> <li>• 1:5 Increasing • For FG burs</li> </ul>	

Viva-mate 3 Motor E015-002  
 • Micromotor 2,000~30,000min<sup>-1</sup>  
 • Includes straight cord

FC-25 Z018  
 • Foot Control (Optional)

### Viva-mate 3 Complete Pack

Viva-mate 3 (120V)	Y103-631
Viva-mate 3 (230V)	Y103-632

- Contents
- Control unit • Micromotor
  - Battery charger
  - Straight handpiece (EX-6B)
  - Carrying case



**Carrying Case**  
Included in the complete pack.

### Control panel



**NAKANISHI INC.**  
www.nsk-inc.com  
700 Shimohinata Kanuma-shi  
Tochigi 322-8666,  
Japan

TEL : +81 (0)289-64-3380  
FAX : +81 (0)289-62-5636

**NSK Europe GmbH**   
www.nsk-europe.de  
Westerbachstraße 58  
D-60489 Frankfurt,  
Germany

TEL : +49 (0)69 74 22 99 0  
FAX : +49 (0)69 74 22 99 29

**NSK United Kingdom Ltd.**  
www.nsk-uk.com  
Office 5, Gateway1000,  
Arlington Business Park,  
Whittle Way, Stevenage, SG1 2FP, UK

TEL : +44 (0) 1438 310670  
FAX : +44 (0) 1438 310671

**NSK Oceania Pty Ltd.**  
www.nsk-inc.com  
Unit 22, 198-222 Young St.  
Waterloo, Sydney,  
NSW 2017 Australia

TEL : +61 (0)2 8306 3000  
FAX : +61 (0)2 8306 3030

**NSK America - Canadian Branch**  
www.nsk-inc.com  
1054 Centre St. - Suite 283  
Thornhill, Ontario, Canada  
L4J 8E5

TEL/FAX : +1 905 780 6757  
Toll Free : +1 866 780 6757

**NSK Middle East**  
www.nsk-inc.com  
Dubai Airport Free Zone,  
Room 2W-310,  
3F West Wing, Dubai, UAE

TEL : +971 (0)4-299-6644  
FAX : +971 (0)4-299-6664