

NSK

ELECTRIC MICROMOTOR

Ti-Max M40

bright, powerful and reliable



- *Clearer*
- *Durable*
- *Safer*



Powerful Partners®

Ti-Max M40

ISO E-Type Internal Spray

The NSK M40 micromotor incorporates advanced coreless DC motor technology. This technology delivers smooth, quiet operation, powerful torque and minimal heat generation despite the M40's compact size.

Advantages of using 



Clearer : LEDs generate natural daylight quality light which assures clearer vision when compared to halogen light.

Durable : LED light sources are longer lasting and more durable than halogen bulbs.

Safer : LED light is safer as it generates less heat than a halogen bulb. Natural daylight quality LED light is easy on the eye.

Features

- The ideal combination of speed and torque
- Minimum heat generation
- Quiet smooth operation from 60 to 40,000 min⁻¹
- LED incorporated in micromotor (M40 LED only)
- Lightweight
- Solid Titanium Body



Optic MODEL M40 LED ORDER CODE E1033 


Optic MODEL M40 ORDER CODE E317

Non-Optic MODEL M40N ORDER CODE E338


Easy to connect to Bien-Air® 4VLM, 4VR 400 and B4VLR hoses



Ti-Max M40



NSK
M40-CD



Bien-Air®
4VLM / 4VR 400 / B4VLR:
(Made by Bien-Air Dental S.A.)

Bien-Air® is registered trademarks of Bien-Air Dental S.A.

NAKANISHI INC. 
www.nsk-inc.com
700 Shimohinata Kanuma-shi
Tochigi 322-8666,
Japan
TEL : +81 (0)289-64-3390
FAX : +81 (0)289-62-5636

NSK Europe GmbH 
www.nsk-europe.de
Elly-Beinhorn-Strasse 8
65760 Eschborn,
Germany
TEL : +49 (0)6196 77606-0
FAX : +49 (0)6196 77606-29

NSK United Kingdom Ltd
www.nsk-uk.com
Office 5, Gateway1000,
Arlington Business Park, Whittle Way,
Stevenage, SG1 2FP, UK
TEL : +44 (0)1438 310670
FAX : +44 (0)1438 310671

NSK Oceania Pty Ltd
www.nsk-inc.com
Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017, Australia
TEL : +61 (0)2 8306 3000
FAX : +61 (0)2 8306 3030

NSK America Corp
www.nsk-inc.com
700 Cooper Court
Schaumburg, IL 60173,
USA
TEL : +1 847 843 7664
FAX : +1 847 843 7622

NSK Middle East
www.nsk-inc.com
Room 6EA-701, 7th Floor, East Wing No.6
Dubai Airport Free Zone,
PO Box 54316 Dubai, UAE
TEL : +971 (0)4-204-5507
FAX : +971 (0)4-204-5508