

# N-CAM

USB 2.0 Intraoral Camera

## A New Design...

NCam features a brand new look and feel with metal housing in titanium color and a slim distal end designed for easy intraoral viewing.

## A New Function...

NCam has a smart 5 step focus ring that allows target focusing from macro to infinity. It also has a new touch activated capture button which minimizes unnecessary movement of the camera when taking pictures.



## A Proprietary Technology...

Like all our cameras, NCam uses our proprietary 7 element glass optical lens to produce unsurpassed image quality in true-to-life color images.



Winus Technology Co., Ltd

From the leader in Precision Optical Design

New Design & Function has all the features you need to perform a quick and accurate intraoral examination

**Image Sample**



- Slim Distal End allows easier access during examination
- Touch sensitive capture button makes it easier to take pictures
- Easy to Use 5 Step Focus captures crystal clear images without distortion



**5 Boxes containing photos taken in the 5 different modes**



1. Smile

2. Arch

3. 2-3 Teeth

4. 1~2 Teeth


5. Macro Mode

**Technical specifications**

Name	N-Cam
Power	DC 5V
Signal System	1/4" SonySuper CCD
Effective Pixels	768 X 494 / 752 X 582
Horizontal Resolution	470 TV Lines
Angle of View	70 degrees
Light Source	4 High Luminance LED
Focusing Type	5 Step Dial
Magnification Rate	1X ~ 35X
Handpiece Length	229mm
Cable Length	3m

**Recommended Computer Requirements**

- Pentium IV 2GHz or faster
- High Speed USB 2.0 port
- At least 512 MB RAM
- High Resolution Monitor(1024 x 768pixels)
- Higher than Windows XP with Service Pack3



Winus Technology Co., Ltd  
From the leader in Precision Optical Design