


# Wcam

USB 2.0 Intraoral Camera

Motorized 5 step focusing button with indication light..



A Proprietary technology from WINUS  
Best of best all our cameras, Wcam uses our proprietary  
7 element glass optical lens to produce unsurpassed  
Image quality in true-to-life color images.

**PRODENT, s.r.o.**

Senická cesta 92; B. Bystrica  
Tel.: 048-4144596 / 0905 480 173  
E-mail: [prodent@prodent.sk](mailto:prodent@prodent.sk)  
[www.prodent.sk](http://www.prodent.sk)

[www.win-us.com](http://www.win-us.com)




Winus Technology Co., Ltd.

From the leader in Precision Optical Design

New Design & Function has all the features you need to perform a quick and accurate intraoral examination

## Camera Functions

**Wcam**



1. LED & Optics
2. Image capture touch button
3. 5 step motor
4. Cable connection part

5 Boxes containing photos taken in the 5 different modes



1. Smile

2. Arch

3. 2-3 Teeth

4. 1~2 Teeth

5. Macro Mode

## Technical specifications

Name	Wcam
Power	5 V DC
Image Sensor	1/4" Sony Exview HAD CCD II 520,000 pixels
Effective Pixel	NTSC 1020(H) X 508(L) / PAL 1020(H) X 596(L) (610K pixels)
Resolution	520 TV lines
Illumination	Four High luminance white chip LED's
focusing type	5 step motorizing by stepping motor
Size & Weight	225mm, 90g
Cable Length	3 meters

## Recommended Computer Requirements

- Pentium IV 2GHz or Faster
- High Speed USB 2.0 Port
- At least 512 MB RAM, 40 GB of Hard Disk
- High Resolution SVGA Monitor (1024x768 pixels or higher)
- More higher version than Windows XP Home or Pro with Service Pack 2.



**WINUS**

Winus Technology Co., Ltd.

From the leader in Precision Optical Design